

Membrane Solutions

Filtration sets for microbiology

BioSet™ Monitor

MS® disposable sterile BioSet™ Monitor are for contaminants monitoring, microbiological testing and sterility testing in liquid samples like food/beverages from raw materials to finished products.

The easy system of low-cost disposable, ready to use sterile filter sets bring significant convenience to your daily activities in the laboratory and enable you to optimize the analysis time for routine testing.



Each unit consisting of a measured filter funnel, base, pad, membrane, removable lid and plug.

Feature:

- MCE gridded membrane with pad with choice of pore size (0.2, 0.45, 0.8 μm)
- Black or white gridded membrane choice for better contrast
- Ready to use: testing time can be significantly reduced
- All-in-one system: filtration unit easily converts to a petri dish, which can be labeled and incubated for culturing
- Heavy-duty polystyrene (PS), translucent, lightweight, and have reinforced rims and inside fluting
- Heat-resistant, chemically inert hydrophobically non stick
- Removable 100ml graduated plastic funnel
- No risk of phthalate contamination
- Easy handling

Application:

- Total colony counts detection, especially for detection of microorganisms contamination in liquids like beverage, wine and so on
- Qualitative and quantitative analysis of microbiological testing of pharmaceuticals, food, beverages, water and other liquids.



Membrane Solutions LLC

Ordering information

Cat No.	Filter color	Pore size (μm)	Filter diameter (mm)	pcs/pk
BM047022W	white/black grid, 100ml	0.22	47mm	50
BM047045W		0.45	47mm	50
BM047080W		0.80	47mm	50
BM047022B	Black/white grid, 100ml	0.22	47mm	50
BM047045B		0.45	47mm	50
BM047080B		0.80	47mm	50

More information pls visit our website www.membrane-solutions.com.

Or mail us at : info@membrane-solutions.com

Membrane Solutions LLC



Membrane Solutions

USA

Toll free: 866-528-4572

Fax: 732-412-4040

Address: 1108 Spring View Lane Plano, TX
75075

China

Toll free: 400-620-1660

Fax: 0086-21- 51687551

Address: 2202, No. 1759 North
Zhongshan Road, Shanghai

For online orders

please write to

order@membrane-solutions.com